# Quick Guide

Beyond the PCR technology,  
BIOFACT promises the progress for your research.

## LYOFACT™ [Lyophilized Type] OneStep RT-PCR Pre-Mix (h)

[Cat. No. BR531-096]

Contents	BR531-096
OneStep RT-PCR Pre-Mix with dye	Lyophilized PCR Tube (0.2 ml) (8 Strip x 12 EA)

### Product Feature

- Preloaded lyophilisates containing Reverse Transcriptase, H-Star Taq polymerase, dNTP's, enhancer and stabilizer
- Simple, and fast setting up procedure for high yield and repeatable results
- Optimized and ready-to-use mixture of all components for a successful qRT-PCR
- M-MLV(Recombinant) : RNase H<sup>-</sup> activity
- Hot start activity : Yes(Chemical-mediated)
- Application : Gene expression analysis, RNA virus research  
Molecular diagnosis

### PCR Mixture & Cycle

RT-PCR Mixture (Reaction vol. : 30 µl)	
Lyophilized OneStep RT-PCR Pre-Mix	1 tube
Primer F (10 pmole/µl)	1 µl
Primer R (10 pmole/µl)	1 µl
Template RNA	< 200 ng
Add DEPC D.W	to 30 µl

Cycle*	
[2-Step cycling protocol]	[3-Step cycling protocol]
50°C 30 min X 1	50°C 30 min X 1
95°C 15 min X 1	95°C 15 min X 1
95°C 20 sec	95°C 20 sec
Anneal & Extension 40 sec } X 25 ~ 40	AT 20-40 sec } X 25 ~ 40
72°C 5 min X 1	72°C 1min/Kb
8°C ∞	72°C 5 min X 1
	8°C ∞

(Template <200 ng)

### Multiple band are observed from the RT-PCR reaction

- The primer sequence is not specific to the desired target.
- The annealing temperature is no optimal.
- The T<sub>m</sub> value of the primers is low
- Number of cycles is too high.
- Genomic DNA contamination-indicated by higher molecular weight product than expected.
- Multiple cDNA products may be obtained from multigene families or alternatively spliced genes.

### Storage

- 4°C for up to 6 months, or -20°C for up to 12 months.
- Store at 4°C in an aluminium coated bag or on another dry place; humidity < 65 % when sealing is opened.

Please contact us, if you have any question and need help.  
T)1670-5695 www.bio-ft.com info@bio-ft.com

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