

Quick Guide

Beyond the PCR technology,
BIOFACT promises the progress for your research.



BioFACT™

RNase Free Water (DEPC-Treated Water)

[Cat. No. PB141-081, PB141-B500]

| Contents | PB141-081 | PB141-B500 |
|---------------------------|--------------|------------|
| BioFACT™ RNase Free Water | 1 mL x 81 ea | 500 mL |

Product description

RNase Free Water is ready to use and requires no preparation, mixing or autoclaving. RNase Free Water has been autoclaved, filtered through a 0.2 µm filter into the final bottle, and autoclaved.

Quality control

RNase Free Water is tested in the following quality control assays before and after final packaging.

Nonspecific endonuclease activity : A sample is incubated with supercoiled plasmid DNA and analyzed by agarose gel electrophoresis.

Exonuclease activity : A sample is incubated with labeled double-stranded DNA, followed by PAGE analysis.

RNase activity : A sample is incubated with labeled RNA and analyzed by PAGE.



Please contact us,
if you have any question and need help.

T)1670-5695

www.bio-ft.com

info@bio-ft.com



2022. 09. 05 (설명서 개정일)

Related Products

| Group | PCR | Purification | Instrument | 전기영동 |
|-------|---|---|--|---|
| 제품사진 | | | | |
| 제품명 | <ul style="list-style-type: none"> - BioFACT™ Taq^{Basic} DNA Polymerase - BioFACT™ H-Star Taq DNA Polymerase - BioFACT™ F-Star Taq DNA Polymerase - BioFACT™ A-Star Taq DNA Polymerase | <ul style="list-style-type: none"> - HiGene™ Plasmid Mini Prep Kit - HiGene™ Gel & PCR Purification System - HiGene™ Genomic DNA Prep Kit (Column Type, Solution Type) - HiGene™ Total RNA Prep Kit | <ul style="list-style-type: none"> - Hush Run™ Thermal Cycler - SimpliAmp™ Thermal Cycler - Mini Heating Dry Bath - BIOZONE Air Sterilizer | <ul style="list-style-type: none"> - gel-rose™ Agarose - EcoDye™ Nucleic Acid Staining Solution - BioFACT™ 1Kb Plus DNA Ladder - BioFACT™ 100bp Plus DNA Ladder |